

This Page Is Inserted by IFW Operations
and is not a part of the Official Record

BEST AVAILABLE IMAGES


Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

IMAGES ARE BEST AVAILABLE COPY.

**As rescanning documents *will not* correct images,
please do not report the images to the
Image Problem Mailbox.**



Biochemicals and Reagents for Life Science Research

ALPHABETICAL
LIST

BIOACTIVE
PEPTIDES

IMMUNO-
CHEMICALS

MOLECULAR
BIOLOGY

RBI,
NEUROSCIENCE,
SIGNAL
TRANSDUCTION

TISSUE
CULTURE

OTHER
PRODUCT
GROUPS/USP

EQUIPMENT,
BOOKS AND
SUPPLIES

DIAGNOSTIC
KITS AND
REAGENTS

PRODUCT
INDEX



SIGMA®

MOLECULAR BIOLOGY PRODUCTS

SEQUENCING

Molecular Biology Products

PRODUCT NUMBER

PRODUCT NUMBER

(Continuation of)
PRIMER SETS, Fluorescent-Labeled for Automated Sequencing

- P 3348** (-40) M13 Forward Primer Set
(5'-GTT TTC CCA GTC ACG ACG-3')
- P 3223** (-29) M13 Reverse Primer Set
(5'-CAG GAA ACA GCT ATG ACC-3')
- P 3473** Lambda gt10 Forward Primer Set
(5'-AGC AAG TTC AGC CTG GTT AAG-3')
- P 6098** Lambda gt10 Reverse Primer Set
(5'-CTT ATG AGT ATT TCT TCC AGG GTA-3')
- P 6723** T3 Phage Promoter Primer Set
(5'-ATT AAC CCT CAC TAA AGG GA-3')
- P 7348** T7 Phage Promoter Primer Set
(5'-TAA TAC GAC TCA CTA TAG GG-3')
- P 7973** SP6 Promoter Primer Set
(5'-ATT TAG GTG ACA CTA TAG-3')
- P 8598** Poly(T)AGC Primer Set
(5'-TTT TTT TTT TTT TTT TTT TTA GC-3')

NUCLEOTIDE SOLUTIONS FOR RESOLVING SEQUENCE AMBIGUITIES

Used for eliminating band compressions caused by the ability of G residues to form secondary structures which are not fully denatured during electrophoresis. Ref.: 1. Mizusawa, S., et al., Nucl. Acids Res., 14, 1319 (1986).
2. Barr, P.J., et al., Biotechniques, 4, 428 (1986).

- D 5163** 2'-Deoxyinosine 1 μ mole 15.50
5'-Triphosphate Shipped in dry ice
(dTTP)
Sodium Salt
5 mM aqueous solution
[95648-77-4] C₁₀H₁₅N₄O₁₃P₃ FW 492.2 (for free acid)
- D 8783** 7-Deaza-2'-deoxyguanosine 0.5 μ mole 54.60
5'-Triphosphate Shipped in dry ice
(-N⁷-dGTP)
Lithium Salt
10 mM aqueous solution
[101515-08-6] C₁₁H₁₇N₄O₁₃P₃ FW 506.2 (for free acid)

DIDEOXYNUCLEOSIDE TRIPHOSPHATE SOLUTIONS 10 mM solutions, pH 7.0

- D 5413** 2',3'-Dideoxyadenosine 0.5 μ mole 22.55
5'-Triphosphate (ddATP) Shipped in dry ice
Lithium Salt
[93939-70-9]
- D 5538** 2',3'-Dideoxycytidine 5'-Triphosphate (ddCTP) 0.5 μ mole 22.55
Lithium Salt Shipped in dry ice
[93939-77-6]
R: 61 S: 45-36/37/39-23
- D 5663** 2',3'-Dideoxyguanosine 0.5 μ mole 22.55
5'-Triphosphate (ddGTP) Shipped in dry ice
Lithium Salt
[93939-69-6]
- D 1789** 2',3'-Dideoxyinosine 0.5 μ mole 56.30
5'-Triphosphate (ddITP) Shipped in dry ice
Lithium Salt
[93858-64-1]
- D 5288** (3'-Deoxythymidine 5'-Triphosphate; ddTTP) (dTTP) 0.5 μ mole 22.55
Lithium Salt Shipped in dry ice
[93939-78-7]

MANUAL AND AUTOMATED SEQUENCING KITS

SEQ-1
DNA SEQUENCING REAGENTS FOR 1 kit
MAXAM-GILBERT (CHEMICAL DEGRADATION) METHODOLOGY
This DNA sequencing kit contains all necessary buffers and modifying reagents used in the Maxam-Gilbert chemical degradation method of DNA sequencing. NO DEA LICENSE IS REQUIRED FOR THIS KIT.
Kit contains sufficient reagents for approx. 100 reactions of sequencing reactions.

D 5154 Dimethyl Sulfate Stop Solution
Used to stop guanosine modification reaction

D 5279 Dimethyl Sulfate
Modification reagent for guanosine

D 5404 Dimethyl Sulfate Buffer
Buffer for guanosine modification reaction

F 4011 Ferric Chloride, 3.0 M
For disposal of hydrazine

H 1764 Hydrazine Stop Solution
Used to stop pyrimidine modification reaction

H 2761 Hydrazine, Anhydrous
Modification reagent for pyrimidines

P 5881 Piperidine
Cleaves DNA strands at modified bases

F 4636 Formic Acid
Modification reagent for purine bases

S 8388 Sodium Acetate, 0.3 M
Used in ethanol precipitation of DNA

S 8513 Sodium Chloride, 5.0 M
Inhibits modification of thymine by hydrazine

S 8263 Sodium Hydroxide, 5.0 M
For disposal of dimethyl sulfate

W 4502 Water, 18 megohm, 0.2 μ m filtered
For reconstitution of buffers and use in modification reactions

Ref.: Maxam, A.M. and Gilbert, W., Meth. Enzymol. 65, 499 (1980).
R: 11-45-46-26/27-28-34-42/43 S: 45-26-36/37/39-23